



ZIMMER BIOMET  
Your progress. Our promise.®




Thoracolumbar Solutions

# Zyston® Curve

## Interbody Spacer System

Surgical Technique Guide



*The Zyston Curve Interbody System is designed to optimize your procedures through simplified insertion and accurate placement.*

## TABLE OF CONTENTS

System Overview	4
Preoperative Planning and Patient Positioning	5
Exposure and Endplate Preparation	6
Trialing and Implant Selection	7
Implantation	8
Implant Removal	10
Implants	11
Instruments	12
Kit Contents	14
Important Information on the Zyston Curve Interbody Spacer System	17


*Zimmer Biomet Spine does not practice medicine. This technique was developed in conjunction with health care professionals. This document is intended for surgeons and is not intended for laypersons. Each surgeon should exercise his or her own independent judgment in the diagnosis and treatment of an individual patient, and this information does not purport to replace the comprehensive training surgeons have received. As with all surgical procedures, the technique used in each case will depend on the surgeon's medical judgment as the best treatment for each patient. Results will vary based on health, weight, activity and other variables. Not all patients are candidates for this product and/or procedure.*

## SYSTEM OVERVIEW

The Zyston Curve Interbody Spacer System implants and supporting instrumentation were designed to help improve the clinical experience of placing TLIF interbody cages in the correct anatomical location.



The system provides a full array of implant options, featuring a bi-directional tapered leading edge, a large autograft cavity and a uniquely designed controlled articulating mechanism to aid in the placement of the implant.

This surgical technique will provide guidance to the approach-related aspects of the TLIF procedure, as well as describe the functionality of the implants and supporting instrumentation.



FEATURES	BENEFITS
Bi-directional tapered leading edge	Self-distracting, aids in implant insertion and distraction
Controlled articulating mechanism	Access multiple insertion angles; allows implant to pivot <i>in situ</i> (up to 55°); facilitates final positioning with minimal passes through the annulotomy window; makes a robust connection to inserter
Large graft cavity	Provides increased volume for autograft packing to help aid in the fusion process
Low profile implant/instrument interface	Allows for increased visualization during implant insertion, particularly in the medial plane
Multiple footprint options	Facilitates a precise anatomical fit
Line-to-line trials	Reduces intraoperative questions regarding final implant size
Unique placement of tantalum markers	Produces clear radiographic visualization during implant insertion

## PREOPERATIVE PLANNING AND PATIENT POSITIONING




### STEP 1

- The patient should be placed prone, in the appropriate position for a posterior approach, and shall be prepared and draped in a manner consistent with surgical facility protocol.
- Utilizing anterior and posterior fluoroscopic imaging and palpation of the patient anatomy, the affected level is identified and marked appropriately for incisions.

**Note:** *The Zyston Curve Spacer System can be implanted using a traditional open approach or implanted using a minimally invasive approach with the AccuVision System.*

Refer to the AccuVision System surgical technique to learn about proper use of the AccuVision System.


## EXPOSURE AND ENDPLATE PREPARATION



### STEP 2

- Upon proper targeting of the affected level(s), a skin incision is made. The soft tissues are dissected and retracted providing the desired visualization of the bony anatomy.
  - The lateral inferior portion of the inferior facet of the superior vertebrae is removed with an osteotome, bur or kerrison. The capsular portion of the ligamentum flavum is exposed and resected. The superior facet of the inferior vertebrae is resected with an osteotome, bur or kerrison.
  - The neural foramen and central spinal canal are decompressed as necessary.
  - The posterolateral portion of the annulus fibrosus is exposed, and an annulus window is created to gain access to the intervertebral space.
- A discectomy is performed.  
**Note:** A posterior discectomy instrument set can be utilized for decompressive and discectomy procedures.
  - The cartilaginous endplates are removed utilizing the paddle scrapers.  
**Note:** Paddle scrapers are available in 1 mm increments from 6 mm to 18 mm. Assemble the modular T-handle to the quick connect fitting of the shaft prior to use.

## TRIALING AND IMPLANT SELECTION



### STEP 3

At the surgeon's discretion, posterior distraction of the vertebral space may be performed.

#### Assembly of Trial Holder

- Insert the threaded inner shaft into the proximal end of the trial holder; turn the inner shaft clockwise to engage the retainer feature of the inner shaft to the trial holder.
- Guide the appropriate size trial head to the distal end of the trial holder ensuring that the concave side of the trial head is in-line with the finger grips of the trial holder.
- Turn the inner shaft clockwise to engage the threads of the trial head until tight.

**Note:** The trials of the Zyston Curve Spacer System have a fixed insertion angle of 20° and are available in 1 mm increments from 6 mm to 18 mm.

- Insert the trial into the annulotomy window and position within the intervertebral space. Confirm positioning with A/P and lateral fluoroscopy.
- Repeat the trial process until the desired amount of distraction is achieved within the intervertebral space.

The height and length of the implant are determined from the final trial.

**Note:** The Zyston Curve trials match the total height of the implant.

**O.R. Tip:** Do not disconnect trials from trial holder.


## IMPLANTATION




MIS variable angle inserter



10°



Variable angle inserter



65°

## STEP 4

It is recommended to pack the anterior portion of the disc space with autograft prior to placement of the Zyston Curve implant.

The Zyston Curve System comes complete with two inserters specific to the needs of the individual procedure.


A minimally invasive (offset) inserter and straight inserter are provided. Both instruments assemble to the implant and function in the same manner.

### Assembly of Implant Inserter

- Insert the threaded inner shaft into the proximal end of the variable angle inserter; turn the inner shaft clockwise to engage the retainer feature of the inner shaft to the inserter.
- Guide the appropriate size Zyston Curve implant to the variable angle inserter, ensuring that the concave side of the implant is in-line with the finger grips of the variable angle inserter.
- Align the titanium insert of the Zyston Curve implant to the inner shaft of the variable angle inserter, and turn the proximal knob clockwise until tight.

**Note:** The Zyston Curve implant can be positioned at any angle between 10° to 65° from the axis of the inserter shaft.





Straight tamp



Footed tamp



Angled tamp



Starter tamp

## STEP 5

- Using the bone graft mold, place the implant into the corresponding slot and fill the graft cavity with autograft.
- Guide the implant into the intervertebral space, verifying placement via fluoroscopic imaging. The angle of insertion can be customized *in situ* to facilitate final positioning.
- **Turn the knob at the proximal end of the inserter counterclockwise until three audible clicks are heard, then move the inserter shaft to the desired angle, turn the knob clockwise until tight, and continue to insert the implant.**
- Repeat the steps until the implant has reached the desired position.
- Verify positioning with A/P and lateral fluoroscopy.
- Remove the inserter by turning the threaded knob at the proximal end of the inserter counterclockwise until free from the implant (approximately four full revolutions).
- Final positioning of the implant can be achieved by using the starter, angled, straight or footed tamps with gentle force.
- Posterior supplemental fixation is performed. See the individual surgical technique manuals for specific instructions.
- Closure is performed per facility aseptic protocols.

### Cleaning


Prior to cleaning and sterilization, remove the inner shafts from the variable angle inserters and trial holders by turning the inner shaft counterclockwise until the retainer feature is free from the instrument and fully remove the inner shaft.

Please refer to the Zimmer Biomet non-sterile instrument IFU for further reprocessing instructions.

## IMPLANT REMOVAL




Implant remover



Slide hammer



Large slide hammer adapter



Small slide hammer adapter

- Locate the threaded portion of the implant.
- Thread the implant remover into the threaded insert.
- Using gentle force, slowly back out implant from the disc space by using the slotted mallet or the slide hammer.


The slide hammer utilizes two types of adapters to connect to the individual instrumentation:

- MIS inserters “small adapter”.
- Straight inserters “large adapter”.
- Implant remover threads into the proximal end of the remover.
- Assemble the appropriate adapter to the slide hammer by threading the adapter onto the distal end of the instrument until tight.

**Note:** The handle can be used to tighten shaft to adapter.

- “Hook” the adapter onto the groove at the proximal portion of the instrument. Pull the slide hammer proximal as necessary to remove the implant.

## IMPLANTS



**Footprint:** 27 mm long × 10 mm wide  
**Available Heights:** 7 mm–16 mm (1 mm increments)  
 17 mm and 18 mm (Available as special order)  
**Lordosis:** 0° and 6°



**Footprint:** 32 mm long × 10 mm wide  
**Available Heights:** 7 mm–16 mm (1 mm increments)  
 17 mm and 18 mm (Available as special order)  
**Lordosis:** 0° and 6°

### Available Graft Volume By Implant Size (Graft volume depicted in cc's):

	LENGTH			
	27 mm		32 mm	
	PARALLEL	LORDOTIC	PARALLEL	LORDOTIC
7 mm	0.55	0.51	0.75	0.70
8 mm	0.63	0.59	0.85	0.81
9 mm	0.71	0.67	0.96	0.91
10 mm	0.79	0.75	1.06	1.02
11 mm	0.86	0.83	1.17	1.12
12 mm	0.94	0.90	1.28	1.23
13 mm	1.02	0.98	1.38	1.34
14 mm	1.10	1.06	1.49	1.44
15 mm	1.18	1.14	1.60	1.55
16 mm	1.26	1.22	1.70	1.66
17 mm	1.34	1.30	1.81	1.76
18 mm	1.41	1.37	1.91	1.87

## INSTRUMENTS



**27 mm Trials** 6 mm–18 mm  
(1 mm increments) PART NUMBERS  
14-533006–14-533018



**32 mm Trials** 6 mm–18 mm  
(1 mm increments) PART NUMBERS  
14-533026–14-533038



**MIS Variable Angle Inserter** PART NUMBER  
14-533072




**MIS Trial Holder** PART NUMBER  
14-533074



**Variable Angle Inserter** PART NUMBER  
14-533021



**Trial Holder** PART NUMBER  
14-533025



**Bone Graft Mold** PART NUMBER  
14-533044



**Implant Remover** PART NUMBER  
14-533019



**Straight Tamp** PART NUMBER  
14-533023




**Angled Tamp** PART NUMBER  
14-533024




**Starter Tamp** PART NUMBER  
14-533046



**Paddle Scrapers** 6 mm–18 mm  
(1 mm increments) PART NUMBERS  
14-533206–14-533218



**Quick Connect T-handle** PART NUMBER  
14-533202



**Slotted Mallet** PART NUMBER  
14-533203



**Large Slide Hammer Adapter** PART NUMBER  
14-533222



**Small Slide Hammer Adapter** PART NUMBER  
14-533221



**Slide Hammer** PART NUMBER  
14-533204



**Footed Tamp** PART NUMBER  
14-533205

## KIT CONTENTS


### Zyston Curve Parallel Implants 14-533105

DESCRIPTION	QTY	PART NUMBER
Curve Spacer, 27 mm × 7 mm, 0°	1	14-533107
Curve Spacer, 27 mm × 8 mm, 0°	2	14-533108
Curve Spacer, 27 mm × 9 mm, 0°	2	14-533109
Curve Spacer, 27 mm × 10 mm, 0°	2	14-533110
Curve Spacer, 27 mm × 11 mm, 0°	2	14-533111
Curve Spacer, 27 mm × 12 mm, 0°	2	14-533112
Curve Spacer, 27 mm × 13 mm, 0°	2	14-533113
Curve Spacer, 27 mm × 14 mm, 0°	2	14-533114
Curve Spacer, 27 mm × 15 mm, 0°	1	14-533115
Curve Spacer, 27 mm × 16 mm, 0°	1	14-533116
Curve Spacer, 27 mm × 17 mm, 0°	0	14-533117*
Curve Spacer, 27 mm × 18 mm, 0°	0	14-533118*
Curve Spacer, 32 mm × 7 mm, 0°	1	14-533147
Curve Spacer, 32 mm × 8 mm, 0°	2	14-533148
Curve Spacer, 32 mm × 9 mm, 0°	2	14-533149
Curve Spacer, 32 mm × 10 mm, 0°	2	14-533150
Curve Spacer, 32 mm × 11 mm, 0°	2	14-533151
Curve Spacer, 32 mm × 12 mm, 0°	2	14-533152
Curve Spacer, 32 mm × 13 mm, 0°	2	14-533153
Curve Spacer, 32 mm × 14 mm, 0°	2	14-533154
Curve Spacer, 32 mm × 15 mm, 0°	1	14-533155
Curve Spacer, 32 mm × 16 mm, 0°	1	14-533156
Curve Spacer, 32 mm × 17 mm, 0°	0	14-533157*
Curve Spacer, 32 mm × 18 mm, 0°	0	14-533158*

### Zyston Curve Lordotic Implants 14-533125


DESCRIPTION	QTY	PART NUMBER
Curve Spacer, 27 mm × 7 mm, 6°	1	14-533127
Curve Spacer, 27 mm × 8 mm, 6°	2	14-533128
Curve Spacer, 27 mm × 9 mm, 6°	2	14-533129
Curve Spacer, 27 mm × 10 mm, 6°	2	14-533130
Curve Spacer, 27 mm × 11 mm, 6°	2	14-533131
Curve Spacer, 27 mm × 12 mm, 6°	2	14-533132
Curve Spacer, 27 mm × 13 mm, 6°	2	14-533133
Curve Spacer, 27 mm × 14 mm, 6°	2	14-533134
Curve Spacer, 27 mm × 15 mm, 6°	1	14-533135
Curve Spacer, 27 mm × 16 mm, 6°	1	14-533136
Curve Spacer, 27 mm × 17 mm, 6°	0	14-533137*
Curve Spacer, 27 mm × 18 mm, 6°	0	14-533138*
Curve Spacer, 32 mm × 7 mm, 6°	1	14-533167
Curve Spacer, 32 mm × 8 mm, 6°	2	14-533168
Curve Spacer, 32 mm × 9 mm, 6°	2	14-533169
Curve Spacer, 32 mm × 10 mm, 6°	2	14-533170
Curve Spacer, 32 mm × 11 mm, 6°	2	14-533171
Curve Spacer, 32 mm × 12 mm, 6°	2	14-533172
Curve Spacer, 32 mm × 13 mm, 6°	2	14-533173
Curve Spacer, 32 mm × 14 mm, 6°	2	14-533174
Curve Spacer, 32 mm × 15 mm, 6°	1	14-533175
Curve Spacer, 32 mm × 16 mm, 6°	1	14-533176
Curve Spacer, 32 mm × 17 mm, 6°	0	14-533177*
Curve Spacer, 32 mm × 18 mm, 6°	0	14-533178*

\*Available by special order.



### Zyston Curve Instruments 14-533000

DESCRIPTION	QTY	PART NUMBER
Trial Head, 27 mm × 6 mm, 70°	1	14-533006
Trial Head, 27 mm × 7 mm, 70°	1	14-533007
Trial Head, 27 mm × 8 mm, 70°	1	14-533008
Trial Head, 27 mm × 9 mm, 70°	1	14-533009
Trial Head, 27 mm × 10 mm, 70°	1	14-533010
Trial Head, 27 mm × 11 mm, 70°	1	14-533011
Trial Head, 27 mm × 12 mm, 70°	1	14-533012
Trial Head, 27 mm × 13 mm, 70°	1	14-533013
Trial Head, 27 mm × 14 mm, 70°	1	14-533014
Trial Head, 27 mm × 15 mm, 70°	1	14-533015
Trial Head, 27 mm × 16 mm, 70°	1	14-533016
Trial Head, 27 mm × 17 mm, 70°	1	14-533017
Trial Head, 27 mm × 18 mm, 70°	1	14-533018
Trial Head, 32 mm × 6 mm, 70°	1	14-533026
Trial Head, 32 mm × 7 mm, 70°	1	14-533027
Trial Head, 32 mm × 8 mm, 70°	1	14-533028
Trial Head, 32 mm × 9 mm, 70°	1	14-533029
Trial Head, 32 mm × 10 mm, 70°	1	14-533030




### Zyston Curve Instruments (Continued) 14-533000

DESCRIPTION	QTY	PART NUMBER
Trial Head, 32 mm × 11 mm, 70°	1	14-533031
Trial Head, 32 mm × 12 mm, 70°	1	14-533032
Trial Head, 32 mm × 13 mm, 70°	1	14-533033
Trial Head, 32 mm × 14 mm, 70°	1	14-533034
Trial Head, 32 mm × 15 mm, 70°	1	14-533035
Trial Head, 32 mm × 16 mm, 70°	1	14-533036
Trial Head, 32 mm × 17 mm, 70°	1	14-533037
Trial Head, 32 mm × 18 mm, 70°	1	14-533038
Curve Implant Remover	1	14-533019
Curve Variable Inserter	1	14-533021
Curve Straight Tamp	1	14-533023
Curve Angled Tamp	1	14-533024
Curve Trial Holder	3	14-533025
Curve Bone Mold	1	14-533044
Curve Starter Tamp	1	14-533046
Curve MIS Variable Inserter	1	14-533072
Curve MIS Trial Holder	1	14-533074



## KIT CONTENTS (Continued)

### Zyston Universal Instruments 14-533200



DESCRIPTION	QTY	PART NUMBER
T-handle	1	14-533202
Slotted Mallet	1	14-533203
Slide Hammer	1	14-533204
Footed Tamp	1	14-533205
Paddle Scraper, 6 mm	1	14-533206
Paddle Scraper, 7 mm	1	14-533207
Paddle Scraper, 8 mm	1	14-533208
Paddle Scraper, 9 mm	1	14-533209
Paddle Scraper, 10 mm	1	14-533210
Paddle Scraper, 11 mm	1	14-533211
Paddle Scraper, 12 mm	1	14-533212
Paddle Scraper, 13 mm	1	14-533213
Paddle Scraper, 14 mm	1	14-533214
Paddle Scraper, 15 mm	1	14-533215
Paddle Scraper, 16 mm	1	14-533216
Paddle Scraper, 17 mm	1	14-533217
Paddle Scraper, 18 mm	1	14-533218
Small Slide Hammer Adapter	1	14-533221
Large Slide Hammer Adapter	1	14-533222



## IMPORTANT INFORMATION ON THE ZYSTON CURVE INTERBODY SYSTEM

### Device Description

The Zyston Curve Interbody Spacer System is intended to be inserted into the intervertebral disc space for intervertebral body fusion.

The Zyston Curve Interbody Spacer is designed to restore height and lordotic angle in the spine. The Zyston Curve Interbody Spacer is available in two lengths (27 mm and 32 mm), two lordotic angles (0 and 6 degrees) and in heights from 7 mm to 18 mm in one millimeter increments.

The device is curved (kidney shaped) to match the anatomical shape of the anterior spine with an endplate sparing design to resist subsidence. The spacer has a tapered leading edge to aid in insertion of the implant, and allows the implant to be self-distracting. The spacer has teeth on the endplate-engaging surfaces to provide stability, resist shear and rotational forces, and to help prevent migration of the spacer within the disc space. The open central cavity allows for placement of autograft material allowing for subsequent bone growth through the interior of the device.

The Zyston Curve Spacer is a lumbar spacer with an internal articulating feature to facilitate implant placement. The design incorporates a threaded titanium insert which has the potential to pivot within the body of the spacer. When connected to the inserter, the insert allows the implant to pivot relative to the shaft, enabling the user to guide the implant into the desired position within the intervertebral space.

The PEEK-OPTIMA® material is radiolucent and permits unobstructed radiographic assessment of the fusion mass. However, due to its radiolucency, the Zyston Curve device has tantalum markers within the body of the spacer to help visualize implant orientation within the spine during surgery and postoperatively.

### Indications for Use

The Zyston Curve Interbody Spacer System is indicated for intervertebral body fusion at one level or two contiguous levels in the lumbar spine from L2 to S1 in patients with degenerative disc disease (DDD) with up to Grade 1 spondylolisthesis at the involved level(s). DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment.

The Zyston Curve Interbody Spacer System is designed for use with autograft to facilitate fusion and is intended for use with supplemental fixation systems cleared for use in the lumbar spine. The Zyston Curve Interbody Spacer System may also be implanted using the AccuVision System to provide the surgeon with a minimally invasive approach for posterior or posterolateral spinal surgery.

### Contraindications

Contraindications include, but are not limited to, infection, systemic, spinal or localized; morbid obesity; signs of local inflammation; fever or leukocytosis; metal sensitivity/allergies to the implant materials; any medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count; grossly distorted anatomy due to congenital abnormalities; rapid joint disease, bone absorption, osteopenia, and/or osteoporosis (osteoporosis is a relative contraindication since this condition may limit the degree of obtainable correction, the amount of mechanical fixation, and/or the quality of the bone graft); any case not needing a bone graft and fusion or where fracture healing is not required; any case requiring the mixing of metals from different components; any patient having inadequate tissue coverage over the operative site or where there is inadequate bone stock, bone quality, or anatomical definition; any case not described in the indications; any patient unwilling to cooperate with the postoperative instructions; any time implant utilization would interfere with anatomical structures or expected physiological performance; prior fusion at the level(s) to be treated.

### Warnings

The surgeon should be aware of the following:

- The correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size, shape and design of the implant. The size and shape of the human spine presents limiting restrictions of the size and strength of implants. No implant can be expected to withstand the unsupported stresses of full weight bearing.
- The surgeon must ensure that all necessary implants and instruments are on hand prior to surgery. The device must be handled and stored carefully, protected from damage, including from corrosive environments. They should be carefully unpacked and inspected for damage prior to use.
- All instruments must be cleaned and sterilized prior to surgery.
- As with all orthopaedic implants, the Zyston Curve Interbody Spacer System should never be reused under any circumstances.
- Proper implant selection and patient compliance to postoperative precautions will greatly affect surgical outcomes. Patients who smoke have been shown to have an increased incidence of nonunion. Therefore, these patients should be advised of this fact and warned of the potential consequences.

## IMPORTANT INFORMATION (Continued)

### Warnings (Continued)

- Postoperative care is important. The patient should be instructed in the limitations of his/her implant and should be cautioned regarding weight bearing and body stress on the appliance prior to secure bone healing.
- The Zyston Curve Interbody Spacer System has not been evaluated for safety and compatibility in the MR environment. The Zyston Curve Interbody Spacer System has not been tested for heating or migration in the MR environment.
- Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without a previous surgery.
- Zimmer Biomet Spine implants should not be used with implants or instruments from another manufacturer for reasons of metallurgy, mechanics and design.
- Based upon dynamic testing results, the physician should consider the levels of implantation, patient weight, patient activity level, other patient conditions, and any other factor which may impact on the performance of the intervertebral body fusion device.

### Precautions

#### Preoperative:

- Only patients that meet the criteria described in the indications should be selected.
- Patient conditions and/or pre-dispositions such as those addressed in the Contraindications section should be avoided.
- Care should be used in the handling and storage of the implant components. The implants should not be scratched or otherwise damaged. Implants and instruments should be protected during storage especially from corrosive environments.
- All instruments should be cleaned and sterilized before use.

#### Intraoperative:

- Any instruction manuals should be carefully followed.
- At all times, extreme caution should be used around the spinal cord and nerve roots. Damage to nerves may occur resulting in a loss of neurological functions.

#### Postoperative:

- The physician's postoperative directions and warnings to the patient and the corresponding patient compliance are extremely important.
- Detailed instructions on the use and limitations of the device should be given to the patient. If partial weight-bearing is recommended or required prior to firm bony union, the patient must be warned that bending, loosening or breakage of the components are complications which can occur as a result of excessive or early weight-bearing

or excessive muscular activity. The risk of bending, loosening, or breakage of an internal fixation device during postoperative rehabilitation may be increased if the patient is active, or if the patient is debilitated, demented, or otherwise unable to use crutches or other such weight supporting devices. The patient should be warned to avoid falls or sudden jolts in spinal position.

- To allow maximum chances for a successful surgical result, the patient or device should not be exposed to mechanical vibrations that may loosen the device construct. The patient should be warned of this possibility and instructed to limit and restrict physical activities, especially lifting, twisting motions and any type of sport participation. The patient should be advised not to smoke or consume alcohol during the bone graft healing process.
- If a nonunion develops or if the components loosen, bend, and/or break, the device(s) should be revised and/or removed immediately before serious injury occurs. Failure to immobilize a delayed or nonunion of bone will result in excessive and repeated stresses on the implant. By the mechanism of fatigue these stresses can cause eventual bending, loosening, or breakage of the device(s). It is important that immobilization of the spinal surgical site be maintained until firm bony union is established and confirmed by radiographic examination. The patient must be adequately warned of these hazards and closely supervised to ensure cooperation until bony union is confirmed.
- Do not reuse implants. While an implant may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant. Do not treat patients with implants/devices that have been even momentarily placed in or used on a different patient.

### Sterilization Recommendations

The Zyston Curve implant is provided sterile. The product is gamma radiation sterilized. The package should be inspected prior to use to ensure the sterile barrier has not been compromised. Do not re-sterilize. Where specified, do not use the device after expiration date.

The Zyston Curve instrumentation is provided non-sterile and must be sterilized prior to use. All packaging materials must be removed prior to use.

Refer to the Zimmer Biomet Non-sterile Instrument IFU for full processing instructions.

Individuals not using the recommended method temperature and time are advised to validate any alternative methods or cycles using an approved method or standard.



**Disclaimer:** This document is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part.



**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician. Rx only. Please refer to the package inserts for important product information, including, but not limited to, indications, contraindications, warnings, precautions, adverse effects, and patient counseling information.

**eLabeling:** The Instructions for Use can be accessed online by visiting the website and using the KEY-CODE provided on the product label and as shown below. Additional languages are also available in electronic format for download. To request a paper copy of the Instructions for Use, contact Zimmer Biomet Spine at the phone number provided.



**Consult Instructions for Use on this website:**

[labeling.zimmerbiomet.com](http://labeling.zimmerbiomet.com)

Key-Code: 14-533099-00



**Manufactured by:**

Zimmer Biomet Spine, Inc.

10225 Westmoor Dr.

Westminster, CO 80021 USA

+1 800.447.3625



**ZIMMER BIOMET**

Your progress. Our promise.®

[800.447.3625/zimmerbiomet.com](http://800.447.3625/zimmerbiomet.com)

©2021 Zimmer Biomet Spine, Inc. All rights reserved.

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet Spine, Inc. or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet Spine. This material is intended for health care professionals, the Zimmer Biomet Spine sales force and authorized representatives. Distribution to any other recipient is prohibited. Check for country product clearances and reference product specific instructions for use.

0751.2-US-en Issue Date: 2021-10-01