

# In Vitro

## Biomechanical Evaluation of a Lateral Lumbar Interbody Fusion Cage with Integrated Lateral Modular Plate Fixation

Lead Investigators

R. DenHaeze,<sup>1</sup> A. Gandhi,<sup>2</sup> C. Ferry,<sup>2</sup> S. Farmer,<sup>2</sup> R. Porter,<sup>3</sup>

Study Presented at:

The European Academy of Neurosurgeons Annual Meeting<sup>†</sup>

**A biomechanical evaluation was performed on a lateral lumbar interbody fusion (LLIF) cage possessing lateral modular plate fixation (MPF) (Timberline MPF Lateral Fusion System – Highridge Medical, Figure 1.)**

The integrated MPF technology, which has yet to be explored in a lateral interbody fusion model, improves vertebral body stabilization during initial cage placement, ensures physiological cage compression, and affords robust bicortical vertebral purchase while preserving the endplates. Given such inherent stabilization capabilities, it is proposed that the LLIF/MPF cage can diminish the need for extensive invasive posterior fixation. Accordingly, it is hypothesized that the LLIF/MPF cage supplemented with minimally disruptive interspinous process fixation (ISPF) can afford comparable or favorable segmental stability to that of a traditional LLIF cage with supplemental bilateral pedicle screw fixation (BPSF)



**Fig. 1** Timberline<sup>®</sup> MPF Lateral Fusion System, 2-Screw (Left) and 4-Screw (Right)

## Study Objective

The objective of this study was to compare the biomechanical stability of the LLIF /MPF construct with supplemental ISPF to that of a traditional LLIF cage with supplemental PSF.

## Study Methods

### Specimen Preparation


Seven human cadaveric spines (L1-L4) were tested. Osseous integrity was confirmed via DEXA scans and radiographs. The L1 and L4 vertebral bodies were potted. Each spine was first tested in an intact state. A lateral discectomy (L2/L3) was then performed, followed by interactive construct instrumentation and testing.

### Testing Protocol:


A 7.5Nm moment was applied in flexion /extension (FE), lateral bending (LB), and axial rotation (AR) using a six degree-of-freedom kinematics system. Segmental range-of-motion (ROM) was tracked using a motion analysis system. Mean ROM reduction and segmental lordosis following fixation was measured

### Testing Sequence:

1. Intact Spine
2. Timberline MPF + ISPF (Alpine XC Adjustable Fusion System) (Figure 2)
3. Lateral Cage BPSF (Figure 3)



**Fig. 2** Timberline® MPF Cage with 4-Hole Plate + Rigid Interspinous Fixation (Alpine XC™ Adjustable Fusion System)



**Fig. 3** Traditional LLIF Cage + Bilateral Pedicle Screw Fixation (Anteroposterior – Left; Lateral – Right) (Anteroposterior – Left; Lateral Right)

## Results




Fig. 4 Flexion/Extension ROM




Fig. 5 Lateral Bending ROM




Fig. 6 Axial Rotation ROM




Fig. 7 Segmental Lordosis (deg.)

## CONCLUSION

The LLIF cage with integrated MPF (4-screw), when supplemented with interspinous process fixation, can provide robust motion reduction in all principle directions; comparing favorable to a traditional LLIF cage augmented with supplemental bilateral PSF.

## References

1. AXIS Neurosurgery, Williamsville, United States,
  2. ZimVie Spine, Broomfield, United States,
  3. Barrow Neurological Institute, Phoenix, United States
- †. DenHaese RP, Gandhi A, Ferry C, Farmer S, Porter R. In-Vitro Biomechanical Evaluation of a Lateral Lumbar Interbody Fusion Cage with Lateral Modular Plate Fixation. European Academy of Neurosurgeons Annual Meeting, Oct 2015, ePoster Presentation.

For more information, visit [highridgemedical.com](https://highridgemedical.com)



**Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician. Rx Only. Please see the product Instructions for Use for a complete listing of the indications, contraindications, precautions, warnings and adverse effects.

Unless otherwise indicated, all trademarks are the property of Highridge. For additional product information, such as risks and contraindications, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Highridge. HM0555 Rev A US-EN 04.2025 ©2025 Highridge Medical, LLC. All rights reserved.