



Will my insurance cover the cost of the device?

Insurance coverage varies depending on your insurance plan. The SpinalPak Stimulator System is generally recognized by Medicare, Medicaid, workers compensation as well as private and public health plans. ZimVie's Patient Advocacy Group is available to assist you and discuss any questions regarding insurance coverage, deductibles and potential out-of-pocket expenses. You may contact our Patient Advocacy Group at 1-888-236-3652.

What is the process for obtaining a SpinalPak Stimulator System?

You can only obtain the SpinalPak Stimulator System with a doctor's prescription.

Who do I contact if I have questions or need to order supplies?

ZimVie Customer Care Representatives are ready to answer your questions, and can be reached at 1-800-526-2579 extension 6000.

For more information, visit [ZimVie.com](https://www.zimvie.com)

ZimVie

1 Gatehall Drive, Suite 303
Parsippany NJ 07054
Phone: +1-800- 526-2579

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to ZimVie or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of ZimVie. This material is intended for health care professionals. Distribution to any other recipient is prohibited. Confidential information intended solely for ZimVie employees, the ZimVie sales force, and authorized representatives. Distribution to any other recipient is prohibited. Complete prescribing information including full indications, contraindications, warnings and precautions associated with the use of these devices may be found online at [ZimVie.com](https://www.zimvie.com). ZV1684 REV A 03/24 ©2024 ZimVie Inc. All rights reserved.



Biomet® SpinalPak® Non-invasive Spine Fusion Stimulator System

PATIENT INFORMATION BROCHURE



Improving daily life




What is the Biomet SpinalPak Non-invasive Spine Fusion Stimulator System?

Following your spine fusion surgery, your doctor may prescribe the SpinalPak Stimulator System which is a proven, safe and effective nonsurgical adjunctive treatment that helps promote the healing of your fusion.

How does the SpinalPak Stimulator System work?

Spinal fusion joins one or more lumbar vertebrae to eliminate motion, increase stability and to try to reduce pain. Following spinal fusion surgery, in the lower (lumbar) spine, the SpinalPak Stimulator System may be prescribed to assist in healing the fusion by sending electrical impulses directly to the spine that mimics your body's natural healing process. It is portable for ambulatory use and easy to operate, so you'll be able to treat your lumbar spinal fusion while going about your daily routine, or while sleeping. Results will vary due to health, weight, activity and other variables. Not all patients are candidates for this medical device. Only a medical professional can determine the treatment appropriate for your specific condition.



Two lightweight electrodes similar in size to a quarter are placed on your back adjacent to the surgical site. The electrodes are easy to apply and are lightweight. The SpinalPak Stimulator System is battery operated with a rechargeable battery pack. Upon connection of the charged battery pack, the SpinalPak Stimulator System is automatically activated and ready to deliver therapeutic treatment.

How long will I have to treat with the device?

The SpinalPak Stimulator System is designed to deliver 270 days of continuous therapeutic treatment. The recommended treatment time is 24 hours per day. Certain risk factors such as smoking, obesity, diabetes, osteoporosis, undergoing a multi-level surgery or a revision spinal surgery may factor into the duration of your treatment.

How do I know the system is working?

You should not feel the SpinalPak Stimulator System as it delivers the treatment signal to the fusion site; however, some patients may feel a slight sensation. When the system is treating, you'll see a battery symbol and check mark on the stimulator's display screen. Please refer to the patient manual for a description of each symbol that may appear on the display screen. If you have any questions, please contact a ZimVie Customer Care Representative at 1-800-526-2579 extension 6000.

Can I use the system while wearing a back brace over my fusion site?

Yes. The SpinalPak Stimulator System's electrodes can be worn comfortably underneath a back brace.

I have a pacemaker. Can I use the system?

The use of a pacemaker or cardioverter must be assessed on an individual basis. ZimVie recommends that you consult with your cardiologist, who can monitor your pacemaker with an electrocardiogram while you are wearing the SpinalPak Stimulator System.

Is the system MRI safe?

MRI scans and procedures should not be performed until the SpinalPak Stimulator System has been completely removed.

Is the system safe to use during pregnancy?

Use of the SpinalPak Stimulator System during pregnancy has not been evaluated; therefore it is not recommended.

