

# Epic™

## Anterior Thoracolumbar Plate

A comprehensive lumbar plating system for robust anterior fixation



**HIGHRIDGE**

# Intuitive Design And Ease Of Use

The Epic Plate provides fixation for use in anterior thoracolumbar fusion procedures. Features include the plate's intuitive design and ease of use, including the single-step cover plate for screw back-out prevention. The Epic System offers straightforward, easy-to-use instrumentation that works to minimize surgical steps



## Epic Design Rationale

- Strong, low-profile plate
- Lumbar and sacral plate options in a wide range of lengths for optimal anatomic fit
- Fixed and variable screw angles
- Unique screw thread design for excellent bone-screw fixation and plate lagging
- 18° conical range of screw insertion angles
- Single-step cover plate for screw back-out prevention
- Simple, intuitive instrumentation



## Implant Dimensions

PLATE	LENGTH
Lumbar	28 mm   30 mm   32 mm   34 mm   36 mm   38 mm 40mm   42 mm   44 mm   46 mm   48 mm   50 mm
Sacral	28 mm   30 mm   32 mm   34 mm   36 mm   38 mm 40 mm   42 mm   44 mm   46 mm*   48 mm*   50 mm*

SCREW	LENGTH
Variable ø6.0mm	20 mm   24 mm   26 mm*   28 mm   30 mm*   32 mm   36 mm
Fixed ø6.0mm	20 mm   24 mm   26 mm*   28 mm   30 mm*   32 mm   36 mm
Rescue ø7.0mm	24 mm   28 mm   32 mm   36 mm

*\*Non-standard, order separately*

For more information, visit [highridgemedical.com](https://highridgemedical.com)

**Distributed by:**  
Highridge Medical  
10225 Westmoor Drive  
Westminster, CO 80021

**Manufactured by:**  
LDR Medical  
Parc d'entreprises du Grand Troyes  
Quartier Europe de l'Ouest  
5, rue de Berlin  
10300 Sainte-Savine France  
+33 (0)3 25 82 32 63



Unless otherwise indicated, all trademarks are the property of Highridge. For additional product information, such as risks and contraindications, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Highridge. HM0501 rev A US-EN 07.2025 ©2025 Highridge Medical, LLC. All rights reserved.