



A comprehensive lumbar plating system for robust anterior fixation.



# Epic™

## Anterior Thoracolumbar Plate



The Epic Plate provides fixation for use in anterior thoracolumbar fusion procedures. Features include the plate's intuitive design and ease of use, including the single-step cover plate for screw back-out prevention. The Epic System offers straightforward, easy-to-use instrumentation that works to minimize surgical steps



## Epic Design Rationale

- Strong, low-profile plate
- Lumbar and sacral plate options in a wide range of lengths for optimal anatomic fit
- Fixed and variable screw angles
- Unique screw thread design for excellent bone-screw fixation and plate lagging
- 18° conical range of screw insertion angles
- Single-step cover plate for screw back-out prevention
- Simple, intuitive instrumentation




## IMPLANT DIMENSIONS

PLATE	LENGTH
Lumbar	28mm   30mm   32mm   34mm   36mm   38mm   40mm   42mm   44mm   46mm   48mm   50mm
Sacral	28mm   30mm   32mm   34mm   36mm   38mm   40mm   42mm   44mm   46mm*   48mm*   50mm*

SCREW	LENGTH
Variable ø6.0mm	20mm   24mm   26mm*   28mm   30mm*   32mm   36mm
Fixed ø6.0mm	20mm   24mm   26mm*   28mm   30mm*   32mm   36mm
Rescue ø7.0mm	24mm   28mm   32mm   36mm

\*Non-standard, order separately



For more information, visit [ZimVie.com](http://ZimVie.com)

**ZimVie Spine**  
10225 Westmoor Drive  
Westminster, CO 80021  
[ZimVie.com](http://ZimVie.com)

**Manufactured by:**  
LDR Medical  
Parc d'entreprises du Grand Troyes  
Quartier Europe de l'Ouest  
5, rue de Berlin  
10300 Sainte-Savine France  
+33 (0)3 25 82 32 63

**Distributed by:**  
ZimVie Spine  
10225 Westmoor Drive  
Westminster, CO 80021

 **ZimVie**  
Restoring Daily Life.®

©2023 ZimVie, Inc. All rights reserved.

All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to ZimVie, Inc. or one of its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of ZimVie. This material is intended for health care professionals and the ZimVie Spine sales force. Distribution to any other recipient is prohibited. PEEK OPTIMA is a trademark of Invivio Biomaterial Solutions. ZV1233-US-en-Rev0523